


XXL Click'aV Plus™ Ligating System

XXL Ligating Clips & Clip Applicators

Safe and minimally invasive solution even for the most challenging cases

Unique size
New possibilities



PROVIDING SPECIALISTS WITH
INNOVATIVE
HEALTHCARE SOLUTIONS

WITH THE CENTRAL FOCUS ON
PATIENTS' SAFETY AND USERS' COMFORT



XXL Click'aV Plus™ Ligating Clips

All the unique features of Click'aV Plus™ – second generation non-absorbable clips are now available for the thick tissues which no other clip on the market can effectively close. Enlarged teeth provide unparalleled transverse and longitudinal stability on the tissue. Due to the reinforced structure clip keeps dimensional stability and remains safely closed on the thickest structures. With Click'aV Plus™ XXL size surgeon can extend the range of economical benefits in procedures performed on extremely thick tissues by replacing expensive endostaplers or harmonic scalpel with the same safe, but cost effective technique.

Availability of both endo and open surgery appliers allows Click'aV Plus XXL™ clips to be used in wide range of procedures. Endoapplier and clip design is optimized to pass 10 mm access port what is great advantage comparing to 12 or 15 mm required for endostaplers.

KEY FEATURES:

- closes structures between 10 and 22 mm through 10 mm trocar
- unlimited hold-on time
- 4 clips per cartridge
- robust clip design
- low interference with CT scanning
- radiolucent in X-rays
- MRI safe
- self-leading lock
- living hinge
- adhesive backing of the cartridge
- ultimate clip stability in the applier and on the tissue
- color-coded cartridges and appliers for easy identification
- compatible with XXL size open and endoappliers

REF	CLIP SIZE	TISSUE BUNDLE SIZE RANGE (MM)	COLOUR	PCS (per cartridge/ cartridges per box)
0301-10XXL04	XXL	10-22	●	4/8



REF	CLIP SIZE	PRODUCT	PCS PER BOX
0301-04XXL20	XXL	Click'aV® Ligating Clip Applier XXL, 20 cm, open surgery	1
0301-04XXLEN	XXL	Click'aV® Clip Applier XXL, endo	1
0301-04XXLEBN	XXL	Click'aV® Clip Applier XXL, endo, bariatric	1

Grena Ltd.
1000 Great West Road,
Brentford, London TW8 9HH,
United Kingdom

**CUSTOMER SERVICE CENTER
NOTTINGHAM OFFICE, UK**
☎ + 44 115 9704 800
☎ + 44 115 9704 800
✉ export@grena.co.uk
www.grena.co.uk



LONDON CHAMBER
COMMERCE AND INDUSTRY